

DIVISION OF INDUSTRY SERVICES
Plumbing Product Review
P.O. Box 2658
Madison, Wisconsin 53701-2658
TTY: Contact Through Relay

Scott Walker, Governor Dave Ross, Secretary

May 20, 2015

DORNOCH MEDICAL SYSTEMS, INC. JEFF O'BRYAN 200 NW PARKWAY RIVERSIDE MO 64150

Re: Description: HEALTH CARE PLUMBING APPLIANCE
Manufacturer: DORNOCH MEDICAL SYSTEMS, INC.
Product Name: TRANSPOSAL ULTRA (trans. ID 2528063)

Model Number(s): UL-EV100, UL-DU100 ULTRA EVAC & ULTRA DUO

Product File No: 20150081

The specifications and/or plans for this plumbing product have been reviewed and determined to be in compliance with chapters SPS 382 through 384, Wisconsin Administrative Code, and Chapters 145 and 160, Wisconsin Statutes.

The Department hereby issues an approval based on the Wisconsin Statutes and the Wisconsin Administrative Code. This approval is valid until the end of May 2020.

This approval supersedes the approval issued on August 18, 2010 under product file number 20100334.

This approval is contingent upon compliance with the following stipulation(s):

- A reduced pressure principle backflow preventer (ASSE 1013 or CAN/CSA B64.4) or air-gap must be installed on the water supply pipes serving this product.
- All reduced principle backflow preventers shall be:
  - 1. installed in accordance with SPS 382.41 (5) Wis. Adm. Code;
  - 2. registered in accordance with SPS 382.20 (13) Wis. Adm. Code; and
  - 3. tested in accordance with SPS 382.21 (8) Wis. Adm. Code.
- The wastewater from these devices shall be discharged to a properly sized sanitary drain, waste and vent (DWV) system and conform to SPS 382.33 (9) (a) 2. and SPS 382.382.10 (1) (f) Wis. Adm. Code.

Existing DWV piping must be large enough to handle the additional wastewater load generated by these devices.

- The wastewater from these devices shall be discharged to a sanitary drain, waste and vent (DWV) system
  comprised of materials that are stable when in contact with the maximum wastewater temperature generated by
  these devices.
- The discharge from this product may discharge to indirect waste piping with an air break or air gap, or be direct connected to the drain waste and vent system.

SBD-10564-E (N.10/97) File Ref: 15008102.DOC

Dornoch Medical Systems Inc. May 20, 2015 Page 2 of 2

Product File No.: 20150081

The department is in no way endorsing this product or any advertising, and is not responsible for any situation which may result from its use.

Sincerely,

Glen W. Schlueter
Environmental Engineer - Plumbing Product Reviewer
Depart of Safety and Professional Services
Division of Industry Services
Bureau of Technical Services
(608) 267-1401 Phone
(608) 267-9723 Fax
glen.schlueter@wi.gov E-mail