



February 23, 2012

PELICAN WATER SYSTEMS
ENVIRO WATER SOLUTIONS
ROBERT PRENTICE
3060 PERFORMANCE CIRCLE, STE. 2
DELAND FL 32724

Re: Description: WATER TREATMENT DEVICE - POU REVERSE OSMOSIS
Manufacturer: PELICAN WATER SYSTEMS
Product Name: PELICAN PRO SIX STAGE RO SYSTEMS
Model Number(s): PELICAN PRO RO
Product File No: 20120036

The specifications and/or plans for this plumbing product have been reviewed and determined to be in compliance with chapters SPS 382 through 384, Wisconsin Administrative Code, and Chapters 145 and 160, Wisconsin Statutes.

The Department hereby issues an approval based on the Wisconsin Statutes and the Wisconsin Administrative Code. This approval is valid until the end of February 2017.

This approval is contingent upon compliance with the following stipulation(s):

- This product has undergone sufficient testing to document the product's ability to reduce only those contaminants and/or substances as specified in this approval letter when the product is installed and maintained in strict accordance with the manufacturer's published instructions.
- If this approved device is modified or additional assertions of function or performance are made, then this approval shall be considered null and void, unless the change is submitted to the department for review and the approval is reaffirmed.
- The system shall be provided with an in-line total dissolved solids (TDS) monitor, or other acceptable means, to warn the user when the system is not performing its functions. Acceptable alternatives to an in-line TDS monitor include:
 1. terminating the discharge of treated water;
 2. sounding an alarm which is connected to acceptable power source;
 3. flashing a light connected to an acceptable power source;
 4. providing the user with an obvious, readily interpretable, indication of the system's ability to perform (e.g. decreasing the flow rate of treated water by 50% or more for systems making mechanical filtration claims;
 5. providing a sampling service by the manufacturer, either directly or through an authorized dealer, a minimum of once every six months;
 6. providing a sampling kit for analysis of TDS or other appropriate contaminants; or
 7. providing a TDS monitor to measure the product water quality.

Whichever means of performance verification is selected, it shall be clearly described in the owner's manual for this device, and approved for use along with the device.

- These devices will only reduce the concentration of cysts/oocysts at water outlets that are served by the devices. Therefore, using point-of-use devices such as these will not protect all routes of potential exposure. Potentially hazardous exposures to cysts/oocysts will remain possible at unprotected outlets.

The presence of cysts/oocysts strongly suggests that other pathogens (e.g. bacteria, virus) may also be present.

If, by way of reputable water analyses, a water supply is known to contain cysts/oocysts, then all the water entering the residence must be treated at the point-of-entry, using an approved water treatment device, to address all potential routes of exposure thereby providing a biologically safe water supply.

- If the treatment components of this device (e.g., replacement cartridge) are replaced with anything other than those originally approved for use with this device, then this approval shall immediately be considered null and void.

Based on testing data submitted to and reviewed by the department, this approval recognizes that this plumbing product will reduce the concentration of contaminants as specified on pages 1 through 3 of this letter.

**HEALTH EFFECTING INORGANIC CONTAMINANT REDUCTION CAPABILITIES
 PRODUCT FILE NUMBER 20120036
 TABLE 1 OF 3**

Production Rate: 55.3 liters per day (lpd) [14.6 gallons per day (gpd)]

Tested Contaminant	Influent Challenge (mg/l) ¹
Arsenic (As ⁺⁵)	0.050 ± 10%
Fluoride (F ⁻¹)	8.0 ± 10%
Lead (Pb ⁺²)	0.15 ± 10%

Other Conditions: the contaminant reduction performance capabilities displayed for Table 1 of 3 were verified by testing conducted in accordance with NSF *International* Standard 58. To qualify for pentavalent arsenic reduction, the device must reduce then influent challenge concentrations such that all effluent concentrations are ≤ 0.010 mg/l. To qualify for Fluoride reduction, the device must reduce the influent challenge concentrations so that all effluent concentrations are ≤ 1.5 mg/l. To qualify for Lead reduction, the device must reduce the influent challenge concentrations so that all effluent concentrations are ≤ 0.010 mg/l.

* = unless otherwise specified
 ≤ = less than or equal to

± = plus or minus
 1 = mg/l are equivalent to parts per million (ppm)

**HEALTH EFFECTING BIOLOGICAL CONTAMINANT REDUCTION CAPABILITIES
 PRODUCT FILE NUMBER 20120036
 TABLE 2 OF 3**

Production Rate: 55.3 liters per day (lpd) [14.6 gallons per day (gpd)]

Capacity: dependent on the type and quantity of particulate matter present in the influent water; the need for maintenance may be indicated by a significant decrease in flow rate.

Tested Contaminant	Influent Challenge (#/ml)
Cysts/Oocysts ¹	≥ 5.0 x 10 ⁴

Other Conditions: the contaminant reduction performance capabilities displayed for Table 2 of 3 were verified by testing conducted in accordance with NSF *International* Standard 53. To qualify for cyst/oocyst reduction, the device must reduce the influent challenge concentrations by ≥ 99.95% at each sample point.

1 = the specific organisms covered under this testing protocol include cryptosporidium parvum, entamoeba histolytica, giardia lamblia and toxoplasma gondii

#/ml = particles per milliliter

≥ = greater than or equal to

**AESTHETIC CONTAMINANT REDUCTION CAPABILITIES
PRODUCT FILE NUMBER 20120036
TABLE 3 OF 3**

Production Rate: 55.3 liters per day (lpd) [14.6 gallons per day (gpd)]

Tested Contaminant	Influent Challenge (mg/l)¹
Total Dissolved Solids (NaCl)	750 ± 40

Other Conditions: the contaminant reduction performance capabilities displayed for Table 3 of 3 were verified by testing conducted in accordance with NSF *International* Standard 58. To qualify for total dissolved solids reduction performance, the device must reduce the influent challenge concentrations by ≥ 75%.

¹ = milligrams per liter (mg/l) are equivalent to parts per million (ppm)
± = plus or minus

≥ = greater than or equal to

This device was tested under controlled laboratory, or field, conditions. The actual performance of this device for a specific end use installation will vary from the tested conditions based on local factors such as water pressure, water temperature and water chemistry.

The department is in no way endorsing this product or any advertising, and is not responsible for any situation which may result from its use.

Sincerely,

Glen W. Schluter
Engineering Consultant
Safety and Buildings Division
Department of Safety and Professional Services
(608) 267-1401 **Phone**
(608) 267-9566 **Fax**
glen.schluter@wi.gov **Email**
7:30AM - 4:30PM CT **Work Hours**

GWS:gws